

Recombinant Human COMMD10

Catalog No: #GP12418



Package Size: #GP12418-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human COMMD10
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 2-202 amino acids of human COMMD10
Target Name	COMM domain containing 10
Other Names	PTD002
Accession No.	Swissprot:Q9Y6G5Gene Accession:BC005179
Uniprot	Q9Y6G5
GeneID	51397;
Storage	-20~-80°C, pH 7.6 PBS

Background

May modulate activity of cullin-RING E3 ubiquitin ligase (CRL) complexes (PubMed:21778237). May down-regulate activation of NF-kappa-B (PubMed:15799966).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only