

# Recombinant Human ABHD14B

Catalog No: #GP12419



Package Size: #GP12419-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ABHD14B
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-210 amino acids of human ABHD14B
Target Name	abhydrolase domain containing 14B
Other Names	CIB; HEL-S-299
Accession No.	Swissprot:Q96IU4Gene Accession:BC007234
Uniprot	Q96IU4
GeneID	84836;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Has hydrolase activity towards p-nitrophenyl butyrate (in vitro). May activate transcription.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only