

Recombinant Human RILPL2

Catalog No: #GP12420



Package Size: #GP12420-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RILPL2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	Rab interacting lysosomal protein-like 2
Other Names	RLP2
Accession No.	Swissprot:Q969X0Gene Accession:BC013042
Uniprot	Q969X0
GeneID	196383;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein that contains a rab-interacting lysosomal protein-like domain. This protein may be involved in regulating lysosome morphology. This protein may also be a target for the Hepatitis C virus and assist in viral replication. Alternate splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only