

Recombinant Human LYPD4

Catalog No: #GP12421



Package Size: #GP12421-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human LYPD4
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	LY6/PLAUR domain containing 4
Other Names	SMR
Accession No.	Swissprot:Q6UWN0Gene Accession:BC034629
Uniprot	Q6UWN0
GeneID	147719;
Storage	-20~-80°C, pH 7.6 PBS

Background

The specific function of LYPD4 is not yet known.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only