

Recombinant Human RAB39B

Catalog No: #GP12422



Package Size: #GP12422-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RAB39B
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	RAB39B, member RAS oncogene family
Other Names	WSMN; MRX72
Accession No.	Swissprot:Q96DA2Gene Accession:BC009714
Uniprot	Q96DA2
GeneID	116442;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the Rab family of proteins. Rab proteins are small GTPases that are involved in vesicular trafficking. Mutations in this gene are associated with X-linked mental retardation.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only