

Recombinant Human SETD3

Catalog No: #GP12424



Package Size: #GP12424-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SETD3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 97-296 amino acids of human SETD3
Target Name	SET domain containing 3
Other Names	C14orf154
Accession No.	Swissprot:Q86TU7Gene Accession:BC009054
Uniprot	Q86TU7
GeneID	84193;
Storage	-20~-80°C, pH 7.6 PBS

Background

Histone methyltransferase that methylates 'Lys-4' and 'Lys-36' of histone H3 (H3K4me and H3K36me). Acts as a transcriptional activator. Plays an important role in the transcriptional regulation of muscle cell differentiation via interaction with MYOD1.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only