

Recombinant Human PLA2G4D

Catalog No: #GP12427



Package Size: #GP12427-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PLA2G4D
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 508-723 amino acids of human PLA2G4D
Target Name	phospholipase A2, group IVD (cytosolic)
Other Names	cPLA2delta
Accession No.	Swissprot:Q86XP0Gene Accession:BC034571
Uniprot	Q86XP0
GeneID	283748;
Storage	-20~-80°C, pH 7.6 PBS

Background

The phospholipase A2 enzyme family, including PLA2G4D, catalyze the hydrolysis of glycerophospholipids at the sn-2 position and then liberate free fatty acids and lysophospholipids (Chiba et al., 2004 [PubMed 14709560]).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only