

Recombinant Human RAB11A

Catalog No: #GP12428



Package Size: #GP12428-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RAB11A
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	RAB11A, member RAS oncogene family
Other Names	YL8
Accession No.	Swissprot:P62491Gene Accession:BC013348
Uniprot	P62491
GeneID	8766;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene belongs to the Rab family of the small GTPase superfamily. It is associated with both constitutive and regulated secretory pathways, and may be involved in protein transport. Two transcript variants encoding different isoforms have been found for this gene.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only