

## Recombinant Human TINAGL1

Catalog No: #GP12430



Package Size: #GP12430-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human TINAGL1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 250-467 amino acids of human TINAGL1
Target Name	tubulointerstitial nephritis antigen-like 1
Other Names	ARG1; LCN7; LIECG3; TINAGRP
Accession No.	Swissprot:Q9GZM7Gene Accession:BC009048
Uniprot	Q9GZM7
GeneID	64129;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is similar in sequence to tubulointerstitial nephritis antigen, a secreted glycoprotein that is recognized by antibodies in some types of immune-related tubulointerstitial nephritis. Three transcript variants encoding different isoforms have been found for this gene.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only