

Recombinant Human PRELID1

Catalog No: #GP12432



Package Size: #GP12432-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PRELID1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 36-219 amino acids of human PRELID1
Target Name	PRELI domain containing 1
Other Names	PX19; PRELI; SBBI12; CGI-106
Accession No.	Swissprot:Q9Y255Gene Accession:BC000007
Uniprot	Q9Y255
GeneID	27166;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the late embryogenesis abundant motif-containing protein family. The encoded protein is localized to mitochondria and may function as a cytoprotectant by regulating cell death and differentiation. Alternative splicing results in multiple transcript variants encoding different isoforms. Several related pseudogenes have been identified.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only