

Recombinant Human ING5

Catalog No: #GP12447



Package Size: #GP12447-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ING5
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to N terminal 226 amino acids of human ING5
Target Name	inhibitor of growth family, member 5
Other Names	p28ING5
Accession No.	Swissprot:Q8WYH8Gene Accession:BC005370
Uniprot	Q8WYH8
GeneID	84289;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is similar to ING1, a tumor suppressor protein that can interact with TP53, inhibit cell growth, and induce apoptosis. This protein contains a PHD-finger, which is a common motif in proteins involved in chromatin remodeling. This protein can bind TP53 and EP300/p300, a component of the histone acetyl transferase complex, suggesting its involvement in TP53-dependent regulatory pathway.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only