

Recombinant Human RNF114

Catalog No: #GP12448



Package Size: #GP12448-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RNF114
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	ring finger protein 114
Other Names	ZNF313; PSORS12
Accession No.	Swissprot:Q9Y508Gene Accession:BC013695
Uniprot	Q9Y508
GeneID	55905;
Storage	-20~-80°C, pH 7.6 PBS

Background

E3 ubiquitin-protein ligase promoting the ubiquitination and degradation of the CDK inhibitor CDKN1A and probably also CDKN1B and CDKN1C. These activities stimulate cell cycle's G1-to-S phase transition and suppress cellular senescence. May play a role in spermatogenesis.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only