

# Recombinant Human RNF141

Catalog No: #GP12451



Package Size: #GP12451-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |   |
|-----------------------|---|
| Product Name          | Recombinant Human RNF141                |
| Brief Description     | Recombinant Protein                     |
| Immunogen Description | Full length fusion protein              |
| Target Name           | ring finger protein 141                 |
| Other Names           | ZFP26; ZNF230                           |
| Accession No.         | Swissprot:Q8WVD5Gene Accession:BC018104 |
| Uniprot               | Q8WVD5                                  |
| GeneID                | 50862;                                  |
| Storage               | -20~-80°C, pH 7.6 PBS                   |

## Background

The protein encoded by this gene contains a RING finger, a motif known to be involved in protein-DNA and protein-protein interactions. Abundant expression of this gene was found in the testicular tissue of fertile men, but was not detected in azoospermic patients. Studies of the mouse counterpart suggest that this gene may function as a testis specific transcription factor during spermatogenesis.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only