

Recombinant Human LYPLA1

Catalog No: #GP12452



Package Size: #GP12452-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human LYPLA1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	lysophospholipase I
Other Names	APT1; LPL1; APT-1; LPL-I; hAPT1
Accession No.	Swissprot:O75608Gene Accession:BC008652
Uniprot	O75608
GeneID	10434;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the alpha/beta hydrolase superfamily. The encoded protein functions as a homodimer, exhibiting both depalmitoylating as well as lysophospholipase activity, and may be involved in Ras localization and signaling. Alternate splicing results in multiple transcript variants. Pseudogenes of this gene have been defined on chromosomes 4, 6, and 7.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only