

Recombinant Human LYPLA2

Catalog No: #GP12453



Package Size: #GP12453-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human LYPLA2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	lysophospholipase II
Other Names	APT2; APT-2; DJ886K2.4
Accession No.	Swissprot:O95372Gene Accession:BC017034
Uniprot	O95372
GeneID	11313;
Storage	-20~-80°C, pH 7.6 PBS

Background

Lysophospholipases are enzymes that act on biological membranes to regulate the multifunctional lysophospholipids. There are alternatively spliced transcript variants described for this gene but the full length nature is not known yet.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only