

Recombinant Human LRRC28

Catalog No: #GP12457



Package Size: #GP12457-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human LRRC28
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-93 amino acids of human LRRC28
Target Name	leucine rich repeat containing 28
Accession No.	Swissprot:Q86X40Gene Accession:BC019704
Uniprot	Q86X40
GeneID	123355;
Storage	-20~-80°C, pH 7.6 PBS

Background

The function of the LRRC28 protein is not known.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only