

Recombinant Human UBE2Z

Catalog No: #GP12464



Package Size: #GP12464-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | Recombinant Human UBE2Z |
| Brief Description | Recombinant Protein |
| Immunogen Description | Fusion protein corresponding to a region derived from 109-354 amino acids of human UBE2Z |
| Target Name | ubiquitin conjugating enzyme E2Z |
| Other Names | USE1; HOYS7 |
| Accession No. | Swissprot:Q9H832Gene Accession:BC015890 |
| Uniprot | Q9H832 |
| GeneID | 65264; |
| Storage | -20~-80°C, pH 7.6 PBS |

Background

This gene encodes an enzyme which ubiquitinates proteins which participate in signaling pathways and apoptosis.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only