

Recombinant Human SPIC

Catalog No: #GP12468



Package Size: #GP12468-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SPIC
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	Spi-C transcription factor (Spi-1/PU.1 related)
Other Names	SPI-C
Accession No.	Swissprot:Q8N5J4Gene Accession:BC032317
Uniprot	Q8N5J4
GeneID	121599;
Storage	-20~-80°C, pH 7.6 PBS

Background

Controls the development of red pulp macrophages required for red blood cells recycling and iron homeostasis. Transcription factor that binds to the PU-box, a purine-rich DNA sequence (5'-GAGGA[AT]-3') that can act as a lymphoid-specific enhancer. Regulates VCAM1 gene expression (By similarity).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only