

# Recombinant Human ASB9

Catalog No: #GP12471



Package Size: #GP12471-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ASB9
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	ankyrin repeat and SOCS box containing 9
Accession No.	Swissprot:Q96DX5Gene Accession:BC001244
Uniprot	Q96DX5
GeneID	140462;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of the ankyrin repeat and suppressor of cytokine signaling (SOCS) box protein family. Members of this family can interact with the elongin B-C adapter complex via their SOCS box domain and further complex with the cullin and ring box proteins to form E3 ubiquitin ligase complexes. They may function to mediate the substrate-recognition of the E3 ubiquitin ligases. A transcribed pseudogene of this gene has been identified on chromosome 15. Multiple transcript variants encoding different isoforms have been found for this gene.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only