

# Recombinant Human HDHD2

Catalog No: #GP12478



Package Size: #GP12478-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human HDHD2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	haloacid dehalogenase-like hydrolase domain containing 2
Other Names	HEL-S-301; 3110052N05Rik
Accession No.	Swissprot:Q9H0R4Gene Accession:BC038100
Uniprot	Q9H0R4
GeneID	84064;
Storage	-20~-80°C, pH 7.6 PBS

## Background

HDHD2 (haloacid dehalogenase-like hydrolase domain containing 2) is also known as DKFZp564D1378 and is a 259 amino acid protein that is expressed as two isoforms produced by alternative splicing. HDHD2 belongs to the HAD-like hydrolase superfamily, which contains a group of hydrolase enzymes

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only