

Recombinant Human ZNF473

Catalog No: #GP12482



Package Size: #GP12482-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ZNF473
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to N terminal 200 amino acids of human ZNF473
Target Name	zinc finger protein 473
Other Names	ZN473; HZFP100
Accession No.	Swissprot:Q8WTR7Gene Accession:BC018612
Uniprot	Q8WTR7
GeneID	25888;
Storage	-20~-80°C, pH 7.6 PBS

Background

Involved in histone 3'-end pre-mRNA processing by associating with U7 snRNP and interacting with SLBP/pre-mRNA complex. Increases histone 3'-end pre-mRNA processing but has no effect on U7 snRNP levels, when overexpressed. Required for cell cycle progression from G1 to S phases.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only