

## Recombinant Human ZFP36

Catalog No: #GP12484



Package Size: #GP12484-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ZFP36
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 127-326 amino acids of human ZFP36
Target Name	ZFP36 ring finger protein
Other Names	TTP; G0S24; GOS24; TIS11; NUP475; zfp-36; RNF162A
Accession No.	Swissprot:P26651Gene Accession:BC009693
Uniprot	P26651
GeneID	7538;
Storage	-20~-80°C, pH 7.6 PBS

## Background

mRNA-binding protein involved in post-transcriptional regulation of AU-rich element (ARE)-containing mRNAs. Acts by specifically binding ARE-containing mRNAs and promoting their degradation. Recruits deadenylase CNOT7 (and probably the CCR4-NOT complex) via association with CNOT1. Plays a key role in the post-transcriptional regulation of tumor necrosis factor (TNF). Plays a key role in the post-transcriptional regulation of tumor necrosis factor (TNF).

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only