

# Recombinant Human ZMAT2

Catalog No: #GP12489



Package Size: #GP12489-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ZMAT2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	zinc finger, matrin-type 2
Other Names	Snu23; Ptg-12; hSNU23
Accession No.	Swissprot:Q96NC0Gene Accession:BC056668
Uniprot	Q96NC0
GeneID	153527;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The function of ZMAT2 remains unknown.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only