

Recombinant Human ZNF19

Catalog No: #GP12491



Package Size: #GP12491-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ZNF19
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to N terminal 250 amino acids of human ZNF19
Target Name	zinc finger protein 19
Other Names	KOX12
Accession No.	Swissprot:P17023Gene Accession:BC047091
Uniprot	P17023
GeneID	7567;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene contains a zinc finger, a nucleic acid-binding domain present in many transcription factors. This gene is located in a region next to ZNF23, a gene also encoding a zinc finger protein, on chromosome 16.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only