

Recombinant Human ZNF213

Catalog No: #GP12494



Package Size: #GP12494-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ZNF213
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 260-459 amino acids of human ZNF213
Target Name	zinc finger protein 213
Other Names	CR53; ZSCAN53; ZKSCAN21
Accession No.	Swissprot:O14771Gene Accession:BC007287
Uniprot	O14771
GeneID	7760;
Storage	-20~-80°C, pH 7.6 PBS

Background

C2H2 zinc finger proteins, such as ZNF213, have bipartite structures in which one domain binds DNA or RNA and the other modulates target gene expression.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only