

Recombinant Human ZNF34

Catalog No: #GP12504



Package Size: #GP12504-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ZNF34
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 22-221 amino acids of human ZNF34
Target Name	zinc finger protein 34
Other Names	KOX32
Accession No.	Swissprot:Q8IZ26Gene Accession:BC004480
Uniprot	Q8IZ26
GeneID	80778;
Storage	-20~-80°C, pH 7.6 PBS

Background

May be involved in transcriptional regulation.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only