

Recombinant Human ZNF43

Catalog No: #GP12510



Package Size: #GP12510-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ZNF43
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 610-809 amino acids of human ZNF43
Target Name	zinc finger protein 43
Other Names	HTF6; KOX27; ZNF39L1
Accession No.	Swissprot:P17038Gene Accession:BC006528
Uniprot	P17038
GeneID	7594;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene belongs to the C2H2-type zinc finger gene family. The zinc finger proteins are involved in gene regulation and development, and are quite conserved throughout evolution. Like this gene product, a third of the zinc finger proteins containing C2H2 fingers also contain the KRAB domain, which has been found to be involved in protein-protein interactions.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only