

Recombinant Human ZNF449

Catalog No: #GP12511



Package Size: #GP12511-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ZNF449
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	zinc finger protein 449
Other Names	ZSCAN19
Accession No.	Swissprot:Q6P9G9Gene Accession:BC065938
Uniprot	Q6P9G9
GeneID	203523;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a nuclear protein that likely functions as a transcription factor. The protein includes an N-terminal SCAN domain, and seven C2H2-type zinc finger motifs.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only