

Recombinant Human ZNF468

Catalog No: #GP12512



Package Size: #GP12512-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ZNF468
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	zinc finger protein 468
Accession No.	Swissprot:Q5VIY5Gene Accession:BC063143
Uniprot	Q5VIY5
GeneID	90333;
Storage	-20~-80°C, pH 7.6 PBS

Background

May be involved in transcriptional regulation.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only