

Recombinant Human ACAT2

Catalog No: #GP12526



Package Size: #GP12526-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ACAT2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 198-397 amino acids of human ACAT2
Target Name	acetyl-CoA acetyltransferase 2
Accession No.	Swissprot:Q9BWD1Gene Accession:BC000408
Uniprot	Q9BWD1
GeneID	39;
Storage	-20~-80°C, pH 7.6 PBS

Background

The product of this gene is an enzyme involved in lipid metabolism, and it encodes cytosolic acetoacetyl-CoA thiolase. This gene shows complementary overlapping with the 3-prime region of the TCP1 gene in both mouse and human. These genes are encoded on opposite strands of DNA, as well as in opposite transcriptional orientation. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only