

Recombinant Human CD14

Catalog No: #GP12527



Package Size: #GP12527-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CD14
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 20-345 amino acids of human CD14
Target Name	CD14 molecule
Accession No.	Swissprot:P08571Gene Accession:BC010507
Uniprot	P08571
GeneID	929;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a surface antigen that is preferentially expressed on monocytes/macrophages. It cooperates with other proteins to mediate the innate immune response to bacterial lipopolysaccharide. Alternative splicing results in multiple transcript variants encoding the same protein.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only