

Recombinant Human TDGF1

Catalog No: #GP12528



Package Size: #GP12528-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TDGF1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 31-150 amino acids of human TDGF1
Target Name	teratocarcinoma-derived growth factor 1
Other Names	CR; CRGF; CRIPTO
Accession No.	Swissprot:P13385Gene Accession:BC022393
Uniprot	P13385
GeneID	6997;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes an epidermal growth factor-related protein that contains a cripto, FRL-1, and cryptic domain. The encoded protein is an extracellular, membrane-bound signaling protein that plays an essential role in embryonic development and tumor growth. Mutations in this gene are associated with forebrain defects. Pseudogenes of this gene are found on chromosomes 2, 3, 6, 8, 19 and X. Alternate splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only