

Recombinant Human YPEL1

Catalog No: #GP12531

Package Size: #GP12531-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human YPEL1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	yippee like 1
Other Names	FKSG3
Accession No.	Swissprot:O60688Gene Accession:BC074501
Uniprot	O60688
GeneID	29799;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene is located in the region associated with DiGeorge syndrome on chromosome 22. The encoded protein localizes to the centrosome and nucleolus and may play a role in the regulation of cell division.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only