

# Recombinant Human YPEL2

Catalog No: #GP12533



Package Size: #GP12533-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human YPEL2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	yippee like 2
Other Names	FKSG4
Accession No.	Swissprot:Q96QA6Gene Accession:BC132884
Uniprot	Q96QA6
GeneID	388403;
Storage	-20~-80°C, pH 7.6 PBS

## Background

YPEL2 (yippee-like 2) belongs to a family of five yippee-like proteins, all of which localize to the centrosome or mitotic spindle and are widely expressed in both adult and fetal tissue. This localization plus the fact that the family of human YPEL proteins share a high degree of sequence homology across species suggests that these proteins may have a conserved function involved in cell division. YPEL2 might be an important factor during the development and malignant transformation of tissues, most notably pancreatic and breast tumors.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only