

# Recombinant Human WDFY4

Catalog No: #GP12542



Package Size: #GP12542-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human WDFY4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 455-654 amino acids of human WDFY4
Target Name	WDFY family member 4
Other Names	C10orf64
Accession No.	Swissprot:Q6ZS81 Gene Accession:BC034937
Uniprot	Q6ZS81
GeneID	57705;
Storage	-20~-80°C, pH 7.6 PBS

## Background

WDFY4 (WDFY family member 4) is a protein-coding gene. Diseases associated with WDFY4 include chediak-higashi syndrome, and systemic lupus erythematosus.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only