

Recombinant Human VGLL2

Catalog No: #GP12553



Package Size: #GP12553-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human VGLL2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	vestigial like family member 2
Other Names	VGL2; VITO1
Accession No.	Swissprot:Q8N8G2Gene Accession:BC069316
Uniprot	Q8N8G2
GeneID	245806;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein with a transcriptional enhancer factor 1 (TEF-1) interaction domain. The encoded protein may act as a co-factor of TEF-1 regulated gene expression during skeletal muscle development. Alternatively spliced transcript variants encoding multiple isoforms have been observed.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only