

Recombinant Human VILL

Catalog No: #GP12554



Package Size: #GP12554-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human VILL
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 657-856 amino acids of human VILL
Target Name	villin-like
Accession No.	Swissprot:O15195Gene Accession:BC004300
Uniprot	O15195
GeneID	50853;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene belongs to the villin/gelsolin family. It contains 6 gelsolin-like repeats and a headpiece domain. It may play a role in actin-bundling.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only