

Recombinant Human UBA5

Catalog No: #GP12566



Package Size: #GP12566-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human UBA5
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 205-404 amino acids of human UBA5
Target Name	ubiquitin-like modifier activating enzyme 5
Other Names	THIFP1; UBE1DC1
Accession No.	Swissprot:Q9GZZ9Gene Accession:BC009737
Uniprot	Q9GZZ9
GeneID	79876;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the E1-like ubiquitin-activating enzyme family. This protein activates ubiquitin-fold modifier 1, a ubiquitin-like post-translational modifier protein, via the formation of a high-energy thioester bond. Two alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. A pseudogene located on chromosome 1 has also been identified.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only