

# Recombinant Human UBASH3A

Catalog No: #GP12569



Package Size: #GP12569-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human UBASH3A
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 424-623 amino acids of human UBASH3A
Target Name	ubiquitin associated and SH3 domain containing A
Other Names	TULA; CLIP4; STS-2; TULA-1
Accession No.	Swissprot:P57075Gene Accession:BC069511
Uniprot	P57075
GeneID	53347;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes one of two family members belonging to the T-cell ubiquitin ligand (TULA) family. Both family members can negatively regulate T-cell signaling. This family member can facilitate growth factor withdrawal-induced apoptosis in T cells, which may occur via its interaction with AIF, an apoptosis-inducing factor. Alternative splicing of this gene results in multiple transcript variants.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only