

Recombinant Human UBE2W

Catalog No: #GP12571



Package Size: #GP12571-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human UBE2W
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	ubiquitin conjugating enzyme E2W (putative)
Other Names	UBC16; UBC-16
Accession No.	Swissprot:Q96B02Gene Accession:BC016326
Uniprot	Q96B02
GeneID	55284;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a nuclear-localized ubiquitin-conjugating enzyme (E2) that, along with ubiquitin-activating (E1) and ligating (E3) enzymes, coordinates the addition of a ubiquitin moiety to existing proteins. The encoded protein promotes the ubiquitination of Fanconi anemia complementation group proteins and may be important in the repair of DNA damage. There is a pseudogene for this gene on chromosome 1. Alternative splicing results in multiple transcript variants.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only