

Recombinant Human UBE2T

Catalog No: #GP12574



Package Size: #GP12574-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human UBE2T
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	ubiquitin conjugating enzyme E2T
Other Names	FANCT; PIG50; HSPC150
Accession No.	Swissprot:Q9NPD8Gene Accession:BC004152
Uniprot	Q9NPD8
GeneID	29089;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene catalyzes the covalent attachment of ubiquitin to protein substrates. Defects in this gene have been associated with Fanconi anemia of complementation group T. Two transcript variants encoding different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only