

Recombinant Human UQCRQ

Catalog No: #GP12591

Package Size: #GP12591-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human UQCRQ
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	ubiquinol-cytochrome c reductase complex III subunit VII
Other Names	QPC; QCR8; QP-C; UQCR7; MC3DN4
Accession No.	Swissprot:O14949Gene Accession:BC001390
Uniprot	O14949
GeneID	27089;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a ubiquinone-binding protein of low molecular mass. This protein is a small core-associated protein and a subunit of ubiquinol-cytochrome c reductase complex III, which is part of the mitochondrial respiratory chain.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only