

Recombinant Human UPK3A

Catalog No: #GP12592



Package Size: #GP12592-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human UPK3A
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 19-287 amino acids of human UPK3A
Target Name	uroplakin 3A
Other Names	UP3A; UPK3; UPIII; UPIIIA
Accession No.	Swissprot:O75631Gene Accession:BC069544
Uniprot	O75631
GeneID	7380;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the uroplakin family, a group of transmembrane proteins that form complexes on the apical surface of the bladder epithelium. Mutations in this gene may be associated with renal adysplasia. Alternatively spliced transcript variants have been described.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only