

Recombinant Human USP6NL

Catalog No: #GP12596



Package Size: #GP12596-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human USP6NL
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 629-828 amino acids of human USP6NL
Target Name	USP6 N-terminal like
Other Names	RNTRE; TRE2NL; USP6NL-IT1
Accession No.	Swissprot:Q92738Gene Accession:BC042943
Uniprot	Q92738
GeneID	9712;
Storage	-20~-80°C, pH 7.6 PBS

Background

Acts as a GTPase-activating protein for RAB5A. Involved in receptor trafficking. In complex with EPS8 inhibits internalization of EGFR.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only