

Recombinant Human UTP11

Catalog No: #GP12602



Package Size: #GP12602-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human UTP11
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	UTP11, small subunit processome component homolog (S. cerevisiae)
Other Names	CGI94; CGI-94; UTP11L
Accession No.	Swissprot:Q9Y3A2Gene Accession:BC005182
Uniprot	Q9Y3A2
GeneID	51118;
Storage	-20~-80°C, pH 7.6 PBS

Background

Involved in nucleolar processing of pre-18S ribosomal RNA.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only