

## Recombinant Human UTP3

Catalog No: #GP12603



Package Size: #GP12603-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human UTP3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 280-479 amino acids of human UTP3
Target Name	UTP3, small subunit processome component homolog (S. cerevisiae)
Other Names	CRL1; CRLZ1; SAS10
Accession No.	Swissprot:Q9NQZ2Gene Accession:BC004546
Uniprot	Q9NQZ2
GeneID	57050;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Essential for gene silencing: has a role in the structure of silenced chromatin. Plays a role in the developing brain (By similarity).

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only