

Recombinant Human UTY

Catalog No: #GP12604



Package Size: #GP12604-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | Recombinant Human UTY |
| Brief Description | Recombinant Protein |
| Immunogen Description | Full length fusion protein |
| Target Name | ubiquitously transcribed tetratricopeptide repeat containing, Y-linked |
| Other Names | UTY1; KDM6AL |
| Accession No. | Swissprot:O14607Gene Accession:BC071744 |
| Uniprot | O14607 |
| GeneID | 7404; |
| Storage | -20~-80°C, pH 7.6 PBS |

Background

This gene encodes a protein containing tetratricopeptide repeats which are thought to be involved in protein-protein interactions. The encoded protein is also a minor histocompatibility antigen which may induce graft rejection of male stem cell grafts. A large number of alternatively spliced transcripts have been observed for this gene, but the full length nature of some of these variants has not been determined.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only