

Recombinant Human ATAD3A

Catalog No: #GP12608



Package Size: #GP12608-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ATAD3A
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 435-634 amino acids of human ATAD3A
Target Name	ATPase family, AAA domain containing 3A
Accession No.	Swissprot:Q9NVI7Gene Accession:BC007803
Uniprot	Q9NVI7
GeneID	55210;
Storage	-20~-80°C, pH 7.6 PBS

Background

ATAD3A and ATAD3B (MIM 612317) are mitochondrial membrane proteins that contribute to the stabilization of large mitochondrial DNA (mtDNA)-protein complexes called nucleoids (He et al., 2007 [PubMed 17210950]).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only