

Recombinant Human ATF7

Catalog No: #GP12610



Package Size: #GP12610-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ATF7
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	activating transcription factor 7
Other Names	ATFA
Accession No.	Swissprot:P17544Gene Accession:BC042363
Uniprot	P17544
GeneID	11016;
Storage	-20~-80°C, pH 7.6 PBS

Background

Plays important functions in early cell signaling. Binds the cAMP response element (CRE) (consensus: 5'-GTGACGT[AG][AG]-3'), a sequence present in many viral and cellular promoters. Activator of the NF-ELAM1/delta-A site of the E-selectin promoter. Has no intrinsic transcriptional activity, but activates transcription on formation of JUN or FOS heterodimers. Also can bind TRE promoter sequences when heterodimerized with members of the JUN family.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only