

Recombinant Human ATL3

Catalog No: #GP12612



Package Size: #GP12612-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ATL3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 342-541 amino acids of human ATL3
Target Name	atlastin GTPase 3
Other Names	HSN1F
Accession No.	Swissprot:Q6DD88Gene Accession:BC077727
Uniprot	Q6DD88
GeneID	25923;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of a family of dynamin-like, integral membrane GTPases. The encoded protein is required for the proper formation of the network of interconnected tubules of the endoplasmic reticulum. Mutations in this gene may be associated with hereditary sensory neuropathy type IF. Alternatively spliced transcript variants that encode distinct isoforms have been described.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only